



# Your PDF Guides


You can read the recommendations in the user guide, the technical guide or the installation guide for SENNHEISER SET 820. You'll find the answers to all your questions on the SENNHEISER SET 820 in the user manual (information, specifications, safety advice, size, accessories, etc.). Detailed instructions for use are in the User's Guide.

**User manual SENNHEISER SET 820**  
**User guide SENNHEISER SET 820**  
**Operating instructions SENNHEISER SET 820**  
**Instructions for use SENNHEISER SET 820**  
**Instruction manual SENNHEISER SET 820**

 **SENNHEISER**

Set 820

Bedienungsanleitung  
Instruction Manual  
Notice d'emploi  
Istruzioni per l'uso  
Instrucciones para el uso  
Gebruiksaanwijzing  
Bruksanvisning



[You're reading an excerpt. Click here to read official SENNHEISER SET 820 user guide](http://yourpdfguides.com/dref/3709644)

<http://yourpdfguides.com/dref/3709644>

*Manual abstract:*

.....  
.....  
.....  
.....

.....  
.....  
.....  
.....

.....  
.....  
.....  
.....

.....  
*...3 English ..*

.....  
.....

.....  
.....  
.....  
.....

.....  
.....  
.....  
.....

.....  
.....  
.....

..... *23 Français* .....

.....  
.....  
.....  
.....

.....  
.....  
.....  
.....

.....  
.....  
.....  
.....

..... *43 Italiano*.....

.....  
.....  
.....

.....  
.....  
.....  
.....

.....  
.....  
.....  
.....



turns on if you pull the earbows apart. After use, the receiver automatically turns off. 28 Adjusting the volume Adjust the receiver's volume with the operating control. If you are not using the stethoset receiver, store it in the transmitter's charging compartment so that the battery can be recharged (see section "Charging the battery").

29 English Changing the transmission channel To avoid interference caused by the simultaneous operation of several transmitters, select a different channel on each transmitter. Several switchable transmission channels are available (CHANNEL). Due to the required approval two channels are selectable: Positions 1 or 3 = channel 1 Position 2 = channel 2 Turning the automatic compression system on Use the COMP switch (compression) at the rear of the transmitter to turn the automatic compression system on. This system reproduces soft passages louder and decreases loud sounds. Use the MO/ST switch to choose between mono and stereo transmission. For mono TVs and when using the MKE 800 TV microphone, the switch must be set to mono (MO). 30 English Adjusting the receiver channel First turn on the transmitter (see section "Where to place the transmitter"). Then put on the stethoset receiver - as soon as the receiver is put on it automatically tunes itself to the first strong signal found (normally the transmitter frequency). If, for some reason, the selected channel is not the correct one (e.g.

@@@The normal charging time is 4-5 hours. A fully charged battery has an operating time of approx. 4 hours. @@@@Delivery includes a pair of replacement ear pads. If necessary, remove the old ear pads and slide on the new ones.

Use a soft cloth for cleaning the units. Do not use abrasive cleaners. 33 WARNING: Volume up? - No! @@@@ATTENTION: Components from other manufacturers (e.g. @@@@Sennheiser batteries can be recycled.

@@@@RR 820 S clip-on receiver Instead of earbows, this receiver features a socket for connecting a headphone or a connecting cable which allows its use with a hearing aid. MKE 800 TV external microphone If your TV doesn't have a headphone output, SCART or RCA socket, you can use your Set 820 together with the MKE 800 TV microphone. Fix the microphone to the loudspeaker of your TV by means of Velcro tape. Connect the microphone connector to the MIC socket on the transmitter and set the transmitter's MO/ST switch to mono (MO). 37 English If problems occur... - Are all jack plugs connected correctly? - Are the receiver and sound source turned on and is the volume of the - receiver and sound source sufficient? - Have you pulled the earbows apart? - Is the battery inserted properly into the receiver? - Have you charged the battery for at least 20 minutes? - Are you too far from the transmitter? (if so, move closer) - Choose a different channel. - Adjust the setting with the TUNE button. - Should the problem persist, please contact your authorized dealer for assistance.

English Never open the unit! 38 Additional information for hearing aid acousticians Adjusting the earbow pressure Loosen the screw at the rear of the stethoset receiver. Move the screw in the long hole to the desired position (continuously adjustable). 1 = Low earbow pressure 2 = High earbow pressure Tighten the screw. 2 1 Adjustments for hearing disabilities Turn the operating control fully counterclockwise (0). Remove the operating control by carefully lifting it over the "guide notch" on the inside. Use the 4 potentiometers on the circuit board to individually adjust the max. volume and the bass response for the left and right headphone system. Bass response L Bass response R Max. volume L Max. volume R 39 English Adjusting the max.

volume: Potentiometer turned fully counterclockwise Potentiometer turned fully clockwise Adjusting the bass response: Potentiometer turned fully counterclockwise Potentiometer turned fully clockwise 0 dB 125 dB attenuation of -25 dB at 100 Hz attenuation 0 dB English The adjustments can also be made when the receiver is being used. @@@@15 kHz 60 dB A rms typ. up to 10 mV 3.5 mm jack, stereo/mono 2.5 mm jack, mono approx.

4 h approx. 10-12 h 12 V DC via NT 820 mains unit approx. @@60 g (without battery) length: 9.0 inch, width: 4.9 inch, height: 1.

0 inch 41 English Manufacturer declarations CE Declaration of Conformity This equipment is in compliance with the essential requirements and other relevant provisions of Directives 1999/5/EC, 89/336/EC or 73/23/EC.



[You're reading an excerpt. Click here to read official SENNHEISER](#)

[SET 820 user guide](#)

<http://yourpdfguides.com/dref/3709644>

*The declaration is available on the internet site at [www.sennheiser.com](http://www.sennheiser.com). Before putting the device into operation, please observe the respective country-specific regulations! English WEEE Declaration Your Sennheiser product was developed and manufactured with highquality materials and components which can be recycled and/or reused. This symbol indicates that electrical and electronic equipment must be disposed of separately from normal waste at the end of its operational lifetime. Please dispose of this product by bringing it to your local collection point or recycling centre for such equipment. This .*



[You're reading an excerpt. Click here to read official SENNHEISER](#)

[SET 820 user guide](#)

<http://yourpdfguides.com/dref/3709644>