



# Your PDF Guides

You can read the recommendations in the user guide, the technical guide or the installation guide for BRAUN PROFESSIONAL CARE 500. You'll find the answers to all your questions on the BRAUN PROFESSIONAL CARE 500 in the user manual (information, specifications, safety advice, size, accessories, etc.). Detailed instructions for use are in the User's Guide.

**User manual BRAUN PROFESSIONAL CARE 500**  
**User guide BRAUN PROFESSIONAL CARE 500**  
**Operating instructions BRAUN PROFESSIONAL CARE 500**  
**Instructions for use BRAUN PROFESSIONAL CARE 500**  
**Instruction manual BRAUN PROFESSIONAL CARE 500**



[You're reading an excerpt. Click here to read official BRAUN PROFESSIONAL CARE 500 user guide](http://yourpdfguides.com/dref/3279369)  
<http://yourpdfguides.com/dref/3279369>

**Manual abstract:**

If the cord is damaged, take the charging unit to an Oral-B Braun Service Centre. A damaged or non-functioning unit should no longer be used. This product is not intended for use by children under age three. This appliance is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by a person responsible for their safety. Do not place the charger in water or other liquid. Do not disassemble the product except when disposing of the battery. When taking out the battery for disposal, use caution not to short the positive (+) and negative (-) terminals. Do not insert any object into any opening of the appliance/charger. Do not touch the power plug with wet hands. When unplugging, always hold the power plug instead of the cord.

Use this product only for its intended use as described in this manual. Do not use attachments which are not recommended by the manufacturer. If you are undergoing treatment for any oral care condition, consult your dental professional prior to use. a b c d e f g Brush head On/off button Handle Charge indicator light Low charge indicator light Basic charger Brush head holder (depending on model) Your toothbrush has a waterproof handle, is electrically safe and designed for use in the bathroom. Plug the basic charger (f) or the assembled charging unit (see picture B, page 2) into an electrical outlet and place the toothbrush handle (c) on the charging unit.

The green charge indicator light (d) flashes while the handle is being charged. Once it is fully charged, the green charge indicator light shines continuously for 5 seconds; then it turns off. A full charge takes at least 17 hours and enables approximately one week of regular brushing (twice a day, 2 minutes). When the rechargeable battery is low, the red low charge indicator light (e) starts flashing after turning off the toothbrush. For everyday use, the toothbrush handle can be stored on the charging unit to maintain it at full power.

overcharging is impossible. Do not press too hard or scrub. In addition you will also hear a different sound while brushing. When pausing longer than 30 seconds the timer resets. In general bleeding should stop after a few days.

We do not recommend using the Braun Oral-B 3D WHITE brush head with braces. Children under 12 years old should not use the Braun Oral-B 3D WHITE brush head. All brush heads feature blue INDICATOR bristles to help you monitor brush head replacement need. With thorough brushing, twice a day for two minutes, the blue color will fade halfway approximately within 3 months, signalling it is time to change the brush head (see picture 1, page 2). After brushing, rinse your brush head or accessory carefully under running water for some seconds with the handle switched on.

Switch off the handle and remove the brush head or accessory. Clean both parts separately under running water then wipe them dry before reassembling the toothbrush. Never place the charger in dishwasher or water; it should be cleaned with a damp cloth only (see picture 6, page 2). Subject to change without notice. In the interest of protecting the environment, please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at an Oral-B Braun Service Centre or at appropriate collection points provided in your country. guarantee We grant 2 years limited guarantee on the product commencing on the date of purchase. The guarantee becomes void if repairs are undertaken by unauthorised persons and if original Braun parts are not used. To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorised Oral-B Braun Service Centre. .



[You're reading an excerpt. Click here to read official BRAUN PROFESSIONAL CARE 500 user guide](http://yourpdfguides.com/dref/3279369)  
<http://yourpdfguides.com/dref/3279369>