# Your PDF Guides

You can read the recommendations in the user guide, the technical guide or the installation guide for SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER. You'll find the answers to all your questions on the SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER in the user manual (information, specifications, safety advice, size, accessories, etc.). Detailed instructions for use are in the User's Guide.

**User manual SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER**  
**User guide SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER**  
**Operating instructions SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER**  
**Instructions for use SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER**  
**Instruction manual SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER**

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**KH 8105 SILVERCREST™**

**i**

Before reading, unfold both pages containing illustrations and familiarise yourself with all functions of the device.

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 [You're reading an excerpt. Click here to read official SILVERCREST KH 8105 FOREHEAD AND EAR THERMOMETER user guide](http://yourpdfguides.com/dref/3035683)  
<http://yourpdfguides.com/dref/3035683>

@@ Press the "EAR" button 3 again to carry out further measurements. 10. The appliance automatically switches off if it is not operated for longer than 1 minute. This application is not stored in the device. In order to return to standard mode, proceed as follows: 1. Press the "FOREHEAD" button 4 and the "EAR" button 3 repeatedly at the same time until the LC display 5 again shows the room temperature, date and time, as well as the symbol (Illustration D). 2. You are now back in standard mode. The integrated watch can be used as a stopwatch in the range of 0. If the sleep mode has been activated, activate the thermometer again by pressing any button before taking a measurement.

2. Press the "FOREHEAD" button 4 and the "EAR" button 3 repeatedly at the same time until the LC display 5 shows the symbol (Illustration F). 3. The stopwatch is started by pressing the "FOREHEAD" button 4 and stopped by pressing it again. You can start a new measurement by pressing the "MEM / SET" button 1 to delete the previously measured time. In order to return to standard mode, proceed as follows: 1. Press the "FOREHEAD" button 4 and the "EAR" button 3 repeatedly at the same time until the LC display 5 again shows the room temperature, date and time, as well as the symbol (Illustration D). 2. You are now in standard mode again. Press and hold down the "EAR" button 3 and the P "FOREHEAD" button 4 at the same time until "°F" appears behind the temperature display.

Switching from "°F" to "°C" extended by the warranty. This also applies to replaced and repaired parts. Any damage or defects present when the device is purchased must be reported immediately upon unpacking, and not later than two days after the date of purchase. All repairs required after the guarantee period has expired are subject to a charge. The packaging is made entirely of recyclable materials, which you may dispose of at local recycling facilities.

@@ the product is subject to Directive 2002 / 96 / EC (WEEE). @@@@ Take used batteries and / or the appliance itself to the collection facilities provided. Batteries should not be disposed of with domestic waste! Take batteries and / or the product itself to the collection facilities provided. In compliance with the following European Directive, the CE symbol 0123 has been applied to the product: Directive for Medical Products 93 / 42 / EG The associated declaration of conformity has been deposited with the importer: KOMPENASS GMBH BURGSTRASSE 21 44867 BOCHUM, GERMANY www. Ta med dessa mätvärden till läkaren, de kan vara till god hjälp för att ställa korrekt diagnos.

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© by ORFGN Marketing